

Dear customer,

WESTAR is the Cyanagen's patented product line dedicated to Western Blotting:

Westar products cover a wide range of possible signal intensity and duration, from the entry level ECL-SUN, to the top of the line, SUPERNOVA.

Each product is at the top of its market segment for performance, and cost to performance ratio.

NEW in our product family: WESTAR HYPERNOVA for proteins with extreme low abundance.

Order-No.	Substrate	Signal intensity	Sensitivity	Protein quantity	Primary Ab dilution range*	Secondary Ab dilution range*	Special price	Special price
XLS063	WESTAR ECL-SUN	Standard	20 pg	High abundance	from 1/100 to 1/5,000	from 1/1,000 to 1/15,000	€ 75,-- for 250 ml	€ 130,-- for 500ml
XLS071	WESTAR Nova 2.0	Medium	2 pg	High abundance	from 1/500 to 1/5,000	from 1/20,000 to 1/100,000	€ 100,-- for 250 ml	€ 165,-- for 500ml
XLS142	WESTAR Antares	High	500 fg	Medium abundance Extended signal duration (20 h)	from 1/1,000 to 1/15,000	from 1/25,000 to 1/150,000	€ 110,-- for 250 ml	€ 200,-- for 500ml
XLS075	WESTAR EtaC Ultra 2.0	Very High	150 fg	Low abundance	from 1/5,000 to 1/50,000	from 1/50,000 to 1/250,000	€ 237,-- for 100ml	€ 399,-- for 200ml
XLS3	WESTAR Supernova	Extreme	50 fg	Very-low abundance	from 1/5,000 to 1/100,000	from 1/100,000 to 1/500,000	€ 284,-- for 100ml	€ 473,-- for 200ml
XL149	WESTAR Hypernova	Very Extreme	>30 fg	extreme-low abundance	few protocol	few protocol	€160,-- for 20ml	

* from 1 mg/ml of antibody stock solution

WESTAR detection reagents are non-isotopic, luminol-based chemiluminescence substrate, designed for the chemiluminescent detection of immobilized proteins and immobilized nucleic acids conjugated with horseradish peroxidase (HRP). Using this method, it is possible to detect membrane immobilized specific antigens, or sequences of nucleic acids, labeled directly with HRP or indirectly with HRP-labeled antibodies/streptavidin.

Use our quote ECL@7BIO to receive our special price until November 30, 2021.

Detailed ordering information and pack sizes can be found on our webpage at www.7bioscience.com

This promotion cannot be combined with any other promotion or discount.

Orders without quote number will not be discounted. Our standard Terms and Conditions do apply.

7Bioscience GmbH
Dekan-Martin-Straße 21
D-79395 Neuenburg
Germany

Tel.: +49 (0) 7631-793 79 80
Fax: +49 (0) 7631-793 79 81
contact@7bioscience.com
www.7bioscience.com

VAT.No: DE298429539
Steuer-Nr. 12179/03828
Company ID: HRB712473
Managing Director: Norbert Gregor

BANK ACCOUNT
Volksbank Freiburg
IBAN: DE09 6809 0000 0037 3632 00
SWIFT-BIC.: GENO DE61